

PRS01

TOGGLE IN TEST PROG.
MD-11-DZPRB-A

EP-DZPRB-A-DL-A
COPYRIGHT © 1977
FICHE 1 OF 1

JUN 1977
digital
MADE IN USA

B01

EOF1DVOVARSEQ

00010000

770608

PDP10 411

IMHOR1DZPRBASEQ

00010000

770608

000000

.REPT 0

1
2
3
4
5
6
7
8
9
10
11
12
13
14
15
16
17
18
19
20
21
22
23
24
25
26
27
28
29
30
31
32
33
34
35
36
37
38
39
40
41
42
43
44
45
46
47
48
49
50
51
52

IDENTIFICATION

PRODUCT CODE:	MAINDEC-11-DZPRB-A-D
PRODUCT NAME:	PR01 TOGGLE IN PROGRAM FOR DL11
DATE:	APRIL 1977
MAINTAINER:	DIAGNOSTIC ENGINEERING
AUTHOR:	JOHN COMEAU

THE INFORMATION IN THIS DOCUMENT IS SUBJECT TO CHANGE WITHOUT NOTICE AND SHOULD NOT BE CONSTRUED AS A COMMITMENT BY DIGITAL EQUIPMENT CORPORATION. DIGITAL EQUIPMENT CORPORATION ASSUMES NO RESPONSIBILITY FOR ANY ERRORS THAT MAY APPEAR IN THIS MANUAL.

THE SOFTWARE DESCRIBED IN THIS DOCUMENT IS FURNISHED TO THE PURCHASER UNDER A LICENSE FOR USE ON A SINGLE COMPUTER SYSTEM AND CAN BE COPIED (WITH INCLUSION OF DIGITAL'S COPYRIGHT NOTICE) ONLY FOR USE IN SUCH SYSTEM, EXCEPT AS MAY OTHERWISE BE PROVIDED IN WRITING BY DIGITAL.

DIGITAL EQUIPMENT CORPORATION ASSUMES NO RESPONSIBILITY FOR THE USE OR RELIABILITY OF ITS SOFTWARE ON EQUIPMENT THAT IS NOT SUPPLIED BY DIGITAL.

COPYRIGHT (C) 1977, BY DIGITAL EQUIPMENT CORPORATION

.ENDR

```

56
57
58      001000
59      177560
60
61
62
63
64
65
66
67
68
69 001000 012703 177560
70 001004 012702 001040
71 001010 005004
72 001012 005012
73 001014 105713
74 001016 100376
75 001020 016300 000002
76 001024 001772
77 001026 105312
78 001030 005204
79 001032 020012
80 001034 001767
81 001036 000000
82      000001
  
```

```

.ABS
=1000
RBUF=177560
:THIS IS THE TOGGLE IN TEST ROUTINE FOR THR PRS01 CONNECTED WITH A DL11.
:IT EXPECTS A TAPE CONSISTING OF NULLS
:FOLLOWED BY A DECREMENTING PATTERN STARTING WITH 377 AND ENDING WITH 1,
:FOLLOWED BY MORE NULLS
:THE ROUTINE WILL HALT AT LOCAION "ERROR" IF IT DETECTS
:A BAD CODE OR A DL11 ERROR BIT
:WHEN THE PROGRAM HALTS, LOCATION 1040 WILL HOLD THE VALUE OF THE
:CHARACTER EXPECTED. R0 WILL HOLD THE VALUE OF THE CHAR ACTUALLY GOTTEN
:BEFORE RUNNING THIS PROGRAM MAKE SURE THAT
:LOCATION 1002 CONTAINS THE DL11 CONTROL/STATUS REGISTER ADDRESS
START:  MOV    #RBUF,R3      :SETUP THE ADDRESS OF THE DL11 BUFFER REGISTER
        MOV    #ERROR+2,R2  :SETUP ADDRESS OF CODE WE EXPECT
        CLR    R4           :ZERO COUNT OF NON-ZERO CHARS
1$:    CLR    (R2)          :SET THE CODE WE SHOULD EXPECT TO 0
3$:    TSTB   (R3)          :GOT A CHAR FROM THE DL11 YET?
        BPL    3$           :IF NOT, GO BACK AND WAIT FOR ONE
        MOV    2(R3),R0     :YES, GOT A CHAR. TAKE IT FROM THE DL11 BUFFER REG
        BEQ    1$           :IF ITS A NULL GO BACK AND GET THE NEXT CHAR
        DECB  (R2)          :NOT A NULL. DECREMENT THE EXPECTED VALUE
        INC   R4           :ADD 1 TO THE COUNT OF NON-ZERO CHARS
        CMP   R0,(R2)       :IS WHAT WE GOT, WHAT WE EXPECTED TO GET?
        BEQ   3$           :IF SO, GO BACK AND KEEP TESTING
ERROR: HALT                :NO. CODE WE GOT AND CODE EXPECTED ARE DIFFERENT
.END
  
```

E01